

varimed®



**PRODUCT CATALOGUE**

**PHLEBOLOGY & LYMPHOLOGY**

**VARİMED MEDİKAL SAN. VE TİC. LTD. ŞTİ.**

# Varimed Classic Anti-Embolism Stockings

Varimed Classic Line Anti Embolism Stockings are designed in order to obtain maximum comfort and benefit to patients. Anti-embolism stocking which are produced as knee-length, thigh length and thigh length with waist belt types, has a measurement system that covers wide measurement range. Also, different colours are given to the heel region of each size of stockings to provide easy separation. Stocking has an inspection hole which controls the blood circulation under the foot fingers.

Thigh length anti-embolism stocking has dot silicone topband which prevents the stocking from slipping and not distorts the compression profile. Also, an elastic tape which doesn't chafe the leg with excellent elasticity and softness is used during the sewing of thigh length with waist belt anti-embolism stocking.

- 18 mmHg, 23 mmHg, 35 mmHg
- White colour
- Colour coded heel
- High wearing comfort and efficacy
- Soft touch
- Extremely durable yarns
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- Perfect sizing. Six sizes and two length
- Skin friendly - Tested to Oeko-Tex Standard 100
- With upper inspection hole and open toe are produced upon request



Compression Class	Stocking Type			Toe Type			Length			Size					
	AD	AG	AG/T	UPPER INS-PECTION HOLE	LOWER INSPECTION HOLE	OPEN	SHORT	NORMAL	LONG	XS	S	M	L	XL	XXL
I (18 mmHg)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
II (23 mmHg)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
III (35 mmHg)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•



LOWER INSPECTION HOLE



UPPER INSPECTION HOLE



OPEN TOE

# Codes And Descriptions Of Varimed Classic Anti Embolism Stockings

## CCL I Light Compression Varimed Classic Anti-Embolism Stockings (18 mmHg)

CODE	PRODUCT
6110	Light Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
6120	Light Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
6130	Light Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
6111	Light Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
6121	Light Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
6131	Light Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
6112	Light Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
6122	Light Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
6132	Light Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking

## CCL II .Medium Compression Varimed Classic Anti-Embolism Stockings (23 mmHg)

CODE	PRODUCT
6210	Medium Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
6220	Medium Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
6230	Medium Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
6211	Medium Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
6221	Medium Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
6231	Medium Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
6212	Medium Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
6222	Medium Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
6232	Medium Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking


## CCL III Strong Compression Varimed Classic Anti-Embolism Stockings (35 mmHg)

CODE	PRODUCT
6310	Strong Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
6320	Strong Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
6330	Strong Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
6311	Strong Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
6321	Strong Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
6331	Strong Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
6312	Strong Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
6322	Strong Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
6332	Strong Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking



## Size Table

Leg Circumference (cm)	Size					
	XS	S	M	L	XL	XXL
cB	18-20	20-23	23-26	26-29	29-32	32-35
cB1	22-27	24-30	27-34	30-37	34-40	38-44
cC	28-34	31-37	34-41	37-44	41-48	44-53
cD	27-33	30-36	33-40	36-43	40-47	43-52
cE	32-38	35-41	38-45	41-48	45-53	49-59
cF	37-47	41-53	47-58	52-62	58-71	64-79
cG	43-57	49-63	53-67	58-73	64-82	72-93
Heel Colour	Beige	Yellow	Green	Red	Blue	Grey



	Length of leg (cm)	
	SHORT	NORMAL
ID	33-37	38-43
lG	62-70	71-82

\*Varimed Anti-Embolism Stockings two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg ( lD is for knee-length stockings, lG is for thigh-length and thigh-length with waist belt stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

\*Anti-Embolism Stockings feature the inspection hole under the toes. Models with inspection hole on the upper side and with open toes are produced upon request.



# Varimed Comfort Anti-Embolism Stockings

Varimed Comfort Line Anti Embolism Stockings are designed in order to obtain maximum comfort and benefit to patients. Anti-embolism stocking which are produced as knee-length, thigh length and thigh length with waist belt types, has a measurement system that covers wide measurement range. Also, different colours are given to the inspection hole of each size of stockings to provide easy separation. Stocking has an inspection hole which controls the blood circulation under the foot fingers. Thigh length anti-embolism stocking has dot silicone topband which prevents the stocking from slipping and not distorts the compression profile. Also, an elastic tape which doesn't chafe the leg with excellent elasticity and softness is used during the sewing of thigh length with waist belt anti-embolism stocking.

- 18 mmHg, 23 mmHg, 35 mmHg
- Ecru colour
- Colour coded inspection hole
- High wearing comfort and efficacy
- Easy put on and take off
- Freedom of movements
- Soft, elastic toe and heel
- Perfect sizing. Six sizes and two length
- Skin friendly – Tested to Oeko-Tex Standard 100
- With upper inspection hole and open toe are produced upon request



Compression Class	Stocking Type			Foot			Length			Size					
	AD	AG	AG/T	Upper Inspection Hole	Lower Inspection Hole	Open	Short	Normal	Long	XS	S	M	L	XL	XXL
I (18 mmHg)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
II (23 mmHg)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
III (35 mmHg)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•



# Product Types and Sizes For Varimed Comfort Anti-Embolism Stockings

## CCL I Light Compression Varimed Classic Anti-Embolism Stockings (18 mmHg)

CODE	PRODUCT
6110E	Light Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
6120E	Light Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
6130E	Light Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
6111E	Light Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
6121E	Light Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
6131E	Light Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
6112E	Light Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
6122E	Light Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
6132E	Light Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking

## CCL II Medium Compression Varimed Classic Anti-Embolism Stockings (23 mmHg)

CODE	PRODUCT
6210E	Medium Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
6220E	Medium Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
6230E	Medium Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
6211E	Medium Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
6221E	Medium Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
6231E	Medium Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
6212E	Medium Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
6222E	Medium Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
6232E	Medium Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking

## CCL III Strong Compression Varimed Classic Anti-Embolism Stockings (35 mmHg)


CODE	PRODUCT
6310E	Strong Compression Knee Length (AD) Upper Inspection Hole Anti-Embolism Stocking
6320E	Strong Compression Thigh Length (AG) Upper Inspection Hole Anti-Embolism Stocking
6330E	Strong Compression Thigh Length with Waist Belt (AG/T) Upper Inspection Hole Anti-Embolism Stocking
6311E	Strong Compression Knee Length (AD) Lower Inspection Hole Anti-Embolism Stocking
6321E	Strong Compression Thigh Length (AG) Lower Inspection Hole Anti-Embolism Stocking
6331E	Strong Compression Thigh Length with Waist Belt (AG/T) Lower Inspection Hole Anti-Embolism Stocking
6312E	Strong Compression Knee Length (AD) Open Toe Anti-Embolism Stocking
6322E	Strong Compression Thigh Length (AG) Open Toe Anti-Embolism Stocking
6332E	Strong Compression Thigh Length with Waist Belt (AG/T) Open Toe Anti-Embolism Stocking



Thigh Length with Waist Belt

## Size Table

Leg Circumference (cm)	Size					
	XS	S	M	L	XL	XXL
cB	18-20	20-23	23-26	26-29	29-32	32-35
cB1	22-27	24-30	27-34	30-37	34-40	38-44
cC	28-34	31-37	34-41	37-44	41-48	44-53
cD	27-33	30-36	33-40	36-43	40-47	43-52
cE	32-38	35-41	38-45	41-48	45-53	49-59
cF	37-47	41-53	47-58	52-62	58-71	64-79
cG	43-57	49-63	53-67	58-73	64-82	72-93
Heel Colour	Beige	Yellow	Green	Red	Blue	Grey



	Length of Leg (cm)	
	SHORT	NORMAL
ID	33-37	38-43
IG	62-70	71-82

\*Varimed Anti-Embolism Stockings two measurements are sufficient to choose the exact size: The ankle circumference (cB) and length of the leg ( lD is for knee-length stockings, lG is for thigh-length and thigh-length with waist belt stockings). The other leg circumference measurements confirm the selected stocking size. Compare the measurements with those in size table. Do not use the stocking in case of patient's whole leg measurements do not correspond to sizes indicated in table.

\*Anti-Embolism Stockings feature the inspection hole under the toes. Models with inspection hole on the upper side and with open toes are produced upon request.



Knee Length



Thigh Length